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Artwork consists of:

- One 4 x 20-inch graphic attached.

Artwork master disk contains the following file(s):

File Name	Description
AW-24892-001_001_01.zip	Source and Supplier print file
AW-24892-001_001_02.pdf	View File


Artwork prints black and Pantone 186 (red)

REV AUTHORED BY G. NIKHIL	DATE JULY/23/2021	HOLOGIC®		SIGNATURES ON FILE
REV DRAFTED BY G. NIKHIL	DATE JULY/23/2021			
PROPRIETARY: This document contains proprietary data of Hologic, Inc. No disclosure, reproduction or use of any part thereof may be made except by written permission from Hologic.		TITLE IFU, THINPREP NUCLEAR STAIN, IVDR, HARD COPY, EN	DOCUMENT NUMBER AW-24892-001	REV 001
REV RELEASE DATE: 4/1/2022		SIZE A	SHEET 1 OF 1	

Before using this document, please consult Agile for the latest revision.

Form ENG-0034-T01, Rev. 006

INTENDED USE

ThinPrep™ Stain Nuclear Stain is intended to be used in a Papanicolaou staining procedure in conjunction with any ThinPrep Stain Rinse Solution, ThinPrep Stain Bluing Solution, ThinPrep Stain Orange G Solution, and ThinPrep Stain EA Solution. For professional use.

SUMMARY AND EXPLANATION

ThinPrep Nuclear Stain is an aqueous solution that serves to stain the nuclei of cells on slide preparations for cytologic evaluation.

PRINCIPLES OF PROCEDURE

When used in conjunction with the ThinPrep staining solutions and cytologic slide preparations, ThinPrep Nuclear Stain will differentially stain the nuclei of cells without interfering with cytoplasmic staining and morphology.

COMPOSITION

An aqueous solution of Acetic acid, Aluminium Sulfate, Ethylene Glycol, Hematoxylin and Sodium Iodate. CAS 107-21-1, CAS 10043-01-3

WARNINGS

Warning. Contains Aluminum Sulfate, Glacial acetic acid, Ethylene glycol. H302. Harmful if swallowed.

For In Vitro Diagnostic Use. Not for external or internal use in humans or animals. May be harmful if swallowed. May cause eye irritation on contact.

PRECAUTIONS

P264. Wash hands thoroughly after handling.
 P280. Wear protective gloves and eye/face protection.
 Avoid contact with eyes. Keep container closed when not in direct use. Use within 60 days after first opening of container. Results with non-ThinPrep cytology samples have not been evaluated. Do not use if primary packaging is damaged.

PRETREATMENT

No reconstitution, mixing or dilution is required.

STORAGE

Store between 15°C (59°F) and 25°C (77°F). Use within 60 days after first opening of container or by the expiration date marked on the container, whichever comes first.

APPEARANCE AND INTEGRITY

Opaque, non-sterile solution.

SPECIMEN COLLECTION AND PREPARATION

For use with gynecologic cytology samples that have been processed on any ThinPrep™ Processor.

PROCEDURE

See the ThinPrep Stain User's Manual for specific staining protocols.

LIMITATIONS OF PROCEDURE

Must be used according to the instructions provided in the ThinPrep Stain User's Manual. ThinPrep Nuclear Stain cannot be substituted with any other solution.

SAFETY AND PERFORMANCE CHARACTERISTICS

When used as directed, ThinPrep Nuclear Stain will stain the nuclei of cells on cytologic slide preparations.

If any serious incident occurs related to this device or any components used with this device, report it to Hologic Technical Support and the competent authority local to the patient and/or user.

Refer to the ThinPrep Imaging System operator's manual.

DISPOSAL

Dispose in accordance with all applicable regulations.

FIRST AID MEASURES

IF SWALLOWED: Call a POISON CENTER or doctor/physician if you feel unwell.

See www.hologicds.com for the entire Safety Data Sheet.

Revision History

Revision	Date	Description
AW-24892-001 Rev. 001	3-2022	Clarify instructions. Add instructions regarding reporting serious incidents.